

4: IDENTIFYING SECURITY DOCUMENTS USING ECCO



Introduction

The most commonly applied method for the forensic examination of paper involves close visual inspection, possibly aided by the wide spectral range of the UV, IR and visible light sources in instruments such as the Foster + Freeman Video Spectral Comparator.

Recent studies have been reported, however, in which paper has been examined using elemental analysis techniques [1, 2, 3]. These studies have tended to make use of the trace elements, e.g. barium (Ba) and strontium (Sr), which occur as impurities in the calcium carbonate (CaCO₃) used as an optical brightener. Security papers, however, tend to have a quite different composition to normal types of paper and are often devoid of optical brighteners. Instead, uncommon elements, such as manganese (Mn) or titanium (Ti), may be present, either having been added intentionally during manufacture or occurring incidentally as a constituent of organic dyes.

In this Application Note, three different types of security paper are examined using the ECCO. Each shows different elemental profiles.

The ECCO relies on the technique of Laser Induced Breakdown Spectroscopy (LIBS) and involves the generation of optical emission spectra from the target material that are characteristic of its elemental composition.

Paper Samples

The study reported here involved subjecting small sections of the inner pages of three different passports to analysis using the ECCO.


Results

ECCO provides the user with a spectral “fingerprint” of the elemental composition of the target material, in this case paper. The heights of the characteristic peaks in the ECCO spectrum are a measure of the concentration of the corresponding elements contained in the material.

The spectra overleaf all show clear differences. In particular, the presence of manganese (Mn) or titanium (Ti) successfully discriminates between the three types of passport.

Passport	Element	
	Manganese	Titanium
German	✓	✗
UK	✗	✓
UK (old style)	✗	✗

Note that the ECCO spectrum of typical A4 office photocopy paper is dominated by peaks associated with the calcium (Ca) in the optical brightener although a small peak due to strontium (Sr) is also visible.



Conclusions

The ECCO provides a fast and effective method for detecting trace elements in paper. The ECCO has been used to discriminate between the different types of papers used in certain passports.

It is likely that the ECCO can also discriminate between other types of security paper, such as that used in banknotes, driving licences and vehicle registration documents.

Unlike other elemental analysis techniques, the ECCO requires very little sample preparation. The instrumentation is cost effective, compact and almost free of maintenance.

References

- [1] L.D. Spence, R. Francis and U. Tinggi, J Forensic Sci. 2002, 47, pages 648 – 651.
- [2] L.D. Spence, A.T. Baker and J.P. Byrne, J Anal At Spectrom 2000, 15, pages 813 – 820.
- [3] A.E. Polk and A.E. Allard, J Forensic Sci., 1977, 22, pages 524 – 533.

Sales Support and Feedback

Foster + Freeman welcome feedback from Customers regarding this product. Please contact one of our offices if you would like to pass on your comments. Foster + Freeman are pleased to offer advice, installation, training and on-site maintenance worldwide for all of their products.

FOSTER + FREEMAN LTD
 25 Swan Lane
 Evesham
 Worcestershire
 WR11 4PE
 UK
 ☎ + 44 (0) 1386 41061
 📠 + 44 (0) 1386 765351
 sales@fosterfreeman.com
 www.fosterfreeman.com

Head Office & UK Sales Office

Foster + Freeman reserve the right to alter the specification of this product, accessories and consumables without prior notice.
 This document contains proprietary information that is protected by copyright.
 All rights are reserved.
 No part of this publication may be reproduced in any form whatsoever without the prior, written permission of Foster + Freeman Ltd.
 Copyright © Foster + Freeman Ltd

FOSTER + FREEMAN
 USA INC
 46030 Manekin Plaza
 Sterling
 Virginia 20166
 USA
 ☎ +1 888 445 5048
 📠 +1 888 445 5049
 usoffice@fosterfreeman.com
 www.fosterfreeman.com

US Sales Office