

3: DISCRIMINATING PRINTER INKS WITH THE FORAM



FORAM Raman Spectral Comparator

The low cost and ready availability of inkjet printers has greatly increased the frequency with which documents produced by these machines are encountered by document examiners. Conventional analytical techniques such as visible/IR absorption which are normally so useful in ink examination are not as effective with printed documents produced by inkjet printers. Other techniques such as chromatography involve the destruction of a small portion of the document.

Whilst the application of Raman and SERRS (Surface Enhanced Resonance Raman Scattering) spectroscopy to the analysis of questioned documents is widely discussed in the scientific literature [1, 2, 3, 4], the application of these techniques to the analysis of black inkjet inks is somewhat limited. Littleford et al [4] have used SERRS spectroscopy to probe the structural changes of the chromophore present in black inkjet inks when deposited onto paper. They also give examples of the types of dye that are likely to be found in inkjet inks.

In this Application Note, we demonstrate the potential of the Foster + Freeman FORAM Raman Spectral Comparator to discriminate black inkjet inks when used in conjunction with the SERRS technique.

Raman spectroscopy involves the scattering of laser light from a target material, the analysis of which provides the user with a spectral "fingerprint" of the molecular composition of the material


Ink samples

The study reported here involved subjecting 15 types of ink to analysis using the FORAM with a laser excitation wavelength of 685 nm. The SERRS technique was effected by applying poly-L-lysine (Sigma-Aldrich) and gold colloid (British Biocell) to the ink mark on the document [5] prior to recording each Raman spectrum. Spectra were baseline-corrected using a propriety fluorescence filter.

Serial	Inkjet Printer	Serial	Inkjet Printer
1	Canon M610	9	HP Deskjet 940c
2	Canon Pixma MP520	10	HP Deskjet 1220c
3	Epson Stylus C64	11	HP Deskjet 6122
4	Epson Stylus C66	12	HP Deskjet F2180
5	Epson Stylus DX7450	13	HP Photosmart 7150
6	HP Business Inkjet 1200	14	Lexmark 4300
7	HP C4180	15	Lexmark Black Inkjet Cartridge #70
8	HP Deskjet 820 Cxi		

Results and Discussion

There are $15 \times 14 / 2 = 105$ sample pairs in the study. Most of the spectral pairs showed clear differences, yielding an overall visual discrimination rate of 84% (88 pairs).



Brunelle and Crawford [6] describe the various types of dyes, solvents, dye complexing agents and surfactants typically found in inkjet inks. The dye component, which is the component expected to give rise to the spectra shown above, is frequently an azo dye with a very broad visible absorption profile. The differences between the dyes are often due to modification or addition of side chain groups [4] to improve properties such as light fastness or solubility. Although it has not been possible to identify the dyes giving rise to the different spectra shown, the small spectral differences observed are consistent with the assertion that the dye molecules have a similar basic molecular structure, but have different side chain groups. Further work is needed to prove this assertion.

Conclusions

The FORAM spectrometer has the ability to discriminate between different types of printer ink when the SERRS technique is applied. A discrimination rate of 84% was achieved. The instrumentation is cost effective, compact and almost free of maintenance.

References

- [1] T.Andermann, Problems of Forensic Sciences, Volume 46 (2001), pp. 335 – 344.
- [2] W.D. Mazzella and P.Buzzini, Forensic Science International 152 (2005), pp. 241-247.
- [3] M.Kunicki, Problems of Forensic Science, Volume 51 (2002), pp. 56 – 70.
- [4] R.E.Littleford, M.P.Hughes, G.Dent, D.Tackley and W.E.Smith, Applied Spectroscopy, Volume 57, Number 8 (2003), pp. 977 – 983.
- [5] E.Wagner and S.Clement, Problems of Forensic Sciences, Volume 46 (2001), Pages 437 – 441.
- [6] “Advances in the forensic analysis and dating of inks”, R.L.Brunelle and K.R.Crawford, published by Charles C Thomas, Springfield, Illinois 2003, pp. 41 – 44.

Sales Support and Feedback

Foster + Freeman welcome feedback from Customers regarding this product. Please contact one of our offices if you would like to pass on your comments. Foster + Freeman are pleased to offer advice, installation, training and on-site maintenance worldwide for all of their products.

FOSTER + FREEMAN LTD
25 Swan Lane
Evesham
Worcestershire
WR11 4PE
UK

☎ + 44 (0) 1386 41061
☎ + 44 (0) 1386 765351
sales@fosterfreeman.com
www.fosterfreeman.com

Head Office & UK Sales Office

Foster + Freeman reserve the right to alter the specification of this product, accessories and consumables without prior notice.

This document contains proprietary information that is protected by copyright.

All rights are reserved.

No part of this publication may be reproduced in any form whatsoever without the prior, written permission of Foster + Freeman Ltd.

Copyright © Foster + Freeman Ltd

FOSTER + FREEMAN
USA INC
46030 Manekin Plaza
Sterling
Virginia 20166
USA

☎ +1 888 445 5048
☎ +1 888 445 5049
usoffice@fosterfreeman.com
www.fosterfreeman.com

US Sales Office